

Anti-PI 3 Kinase p85 beta (8D9) Mouse antibody



Catalog Number: 168120

Product name

Anti-PI 3 Kinase p85 beta (8D9) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

PI 3 Kinase p85 beta (8D9) Mouse monoclonal antibody

Dilutions

WB: 1/1000


Preparation

Antigen: Purified recombinant human PI3 Kinase p85 beta protein fragments expressed in E.coli.

Validations

Formulation

PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of PI3 Kinase p85 beta in C6,K562,MCF7,U87MG,HCT116,3T3 and HeLa cell lysates using PI3 Kinase p85 beta mouse mAb (1:1000 diluted).Predicted band size:85KDa.Observed band size:85KDa.

Clonality

Monoclonal