

Anti-Rb1 Mouse antibody

Catalog Number: 169120

Product name

Anti-Rb1 Mouse antibody

Specificity

Human

Antibody description

Rb1 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human RB1 (AA: 2106-2784) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide and 0.5% protein stabilizer

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

WB, FC


Validations



Flow cytometric analysis of Jurkat cells using RB1 mouse mAb (green) and negative control (red).



Western blot analysis using RB1 mAb against



human RB1 (AA)

Western blot analysis using RB1 mouse mAb against Jurkat (1) and A431 (2) cell lysate.