

# FANCI antibody



Catalog Number: 110520

## Product name

FANCI antibody

## Immunogen

[Human FANCI Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

FANCI Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive IP detected in Jurkat cells. Positive WB detected in HepG2 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 150 kDa

## Preparation

This antibody was obtained by immunization of FANCI recombinant protein (Accession Number: BC004277). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


IP: 1:500-1:5000

IF: 1:10-1:100

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110520(FANCI antibody) at dilution of 1:500




IP Result of anti-FANCI (IP:Catalog No:110520, 4ug; Detection:Catalog No:110520 1:1000) with Jurkat cells lysate 2400ug.

# FANCI antibody



Catalog Number: 110520

---



Immunofluorescent analysis of HeLa cells, using FANCI antibody Catalog No:110520 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).