

FANCI antibody

Catalog Number: 110520

Product name

FANCI antibody

Immunogen

[Human FANCI Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

FANCI Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IP detected in Jurkat cells. Positive WB detected in HepG2 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 150 kDa

Preparation

This antibody was obtained by immunization of FANCI recombinant protein (Accession Number: BC004277). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


IP: 1:500-1:5000

IF: 1:10-1:100

Validations




HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110520(FANCI antibody) at dilution of 1:500



IP Result of anti-FANCI (IP:Catalog No:110520, 4ug; Detection:Catalog No:110520 1:1000) with Jurkat cells lysate 2400ug.

FANCI antibody

Catalog Number: 110520



Immunofluorescent analysis of HeLa cells, using FANCI antibody Catalog No:110520 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).