

# OVCA1 antibody

Catalog Number: 113520

## Product name

OVCA1 antibody

IHC: 1:20-1:200

IF: 1:20-1:200


## Specificity

Human; other species not tested.

## Validations

## Antibody description

OVCA1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IF detected in HeLa cells. Positive IHC detected in human heart tissue. Observed molecular weight by Western-blot: 40-60kd




## Preparation

This antibody was obtained by immunization of OVCA1 recombinant protein (Accession Number: NM\_001383). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113520(OVCA1 antibody) at dilution of 1:300

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human heart using Catalog No:113520(OVCA1 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG



## Applications

ELISA, WB, IHC, IF

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000


Immunohistochemical of paraffin-embedded human heart using Catalog No:113520(OVCA1

# OVCA1 antibody

Catalog Number: 113520

---

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:113520(OVCA1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)