

# RADIL antibody



Catalog Number: 114520

## Product name

RADIL antibody

IP: 1:200-1:2000

IF: 1:10-1:100

## Specificity

Human; other species not tested.

## Antibody description

RADIL Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive IP detected in HepG2 cells. Positive WB detected in HEK-293 cells, HeLa cells, HepG2 cells, human brain tissue, human testis tissue. Observed molecular weight by Western-blot: 117kd

## Preparation

This antibody was obtained by immunization of RADIL recombinant protein (Accession Number: BC004919). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations




HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114520(RADIL antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using RADIL antibody Catalog No:114520 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent

DNA dye).



# RADIL antibody



Catalog Number: 114520

---

IP Result of anti-RADIL (IP:Catalog No:114520,

5ug; Detection:Catalog No:114520 1:800) with  
HepG2 cells lysate 1720ug.