

Anti-ATP1B1 Rabbit antibody



Catalog Number: 167520

Product name

Anti-ATP1B1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ATP1B1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human ATP1B1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions


WB: 1/1000-1/5000

IHC: 1/20-1/50

Validations



Western blot detection of ATP1B1 in Rat Brain lysates using ATP1B1 Rabbit pAb(1:1000 diluted). Predicted band size:35kDa.Observed band size:45-55kDa.



Immunohistochemistry of ATP1B1 in paraffin-embedded Human colon cancer tissue using ATP1B1 Rabbit pAb at dilution 1/20