

AXUD1 antibody



Catalog Number: 108520

Product name

AXUD1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AXUD1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, mouse thymus tissue. Observed molecular weight by Western-blot: 64-70 kDa

Preparation

This antibody was obtained by immunization of AXUD1 recombinant protein (Accession Number: XM_047448724). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117081(AXUD1



antibody) at dilution of 1:500

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117081(AXUD1 antibody) at dilution of 1:500



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:117081(AXUD1 antibody) at dilution of 1:500