

E.coli uidA/gurA/gusA/b1617/JW1609 (His tag) recombinant protein



Catalog Number: 519520

General Information

Gene Name Synonym

Beta-glucuronidase,uidA,GUS

Protein Construction

A DNA sequence encoding the E.coli uidA/gurA/gusA/b1617/JW1609 Met1-Gln603 is expressed with a 6His tag at the N-terminus.

Organism

E.coli

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of PBS,pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



E.coli uidA/gurA/gusA/b1617/JW1609 (His tag)
recombinant protein