

# GDAP1L1 antibody



Catalog Number: 110920

## Product name

GDAP1L1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

GDAP1L1 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells. Positive IHC detected in human gliomas. Observed molecular weight by Western-blot: 42kd

## Preparation

This antibody was obtained by immunization of GDAP1L1 recombinant protein (Accession Number: NM\_024034). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:110920(GDAP1L1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110920(GDAP1L1 antibody) at dilution of 1:50 (under 10x lens)