

# PKNOX1 antibody



Catalog Number: 113920

## Product name

PKNOX1 antibody

IHC: 1:20-1:200

IF: 1:50-1:500

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PKNOX1 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive IF detected in HepG2 cells. Positive WB detected in COLO 320 cells, HepG2 cells. Observed molecular weight by Western-blot: 50-55 kDa

## Preparation

This antibody was obtained by immunization of PKNOX1 recombinant protein (Accession Number: XM\_047440826). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC, WB, IF

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000


## Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:113920(PKNOX1 antibody) at dilution of 1:100 (under 25x lens)



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:113920(PKNOX1 Antibody) at dilution of 1:300



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog

# PKNOX1 antibody



Catalog Number: 113920

---

No:113920(PKNOX1 Antibody) at dilution of 1:50

and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)