

Anti-Phospho-Tau (Ser235) Rabbit antibody



Catalog Number: 166920

Product name

Polyclonal

Anti-Phospho-Tau (Ser235) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB

Antibody description

Phospho-Tau (Ser235) Rabbit polyclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser235 of human Tau

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Western blot detection of Tau (Phospho-Ser235) in Human Hippocampus lysates using Tau (Phospho-Ser235) Rabbit pAb(1:1000 diluted). Predicted band size:79KDa.Observed band size:50-80KDa.

