

# AIRE antibody



Catalog Number: 107920

## Product name

AIRE antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

AIRE Rabbit Polyclonal antibody. Positive IHC detected in human stomach tissue. Positive WB detected in L02 cells, human plasma tissue, human spleen tissue, mouse thymus tissue, U-937 cells. Observed molecular weight by Western-blot: 65-70 kDa

## Preparation

This antibody was obtained by immunization of AIRE recombinant protein (Accession Number: NM\_000658). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107920(AIRE Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human stomach slide using Catalog No:107920(AIRE Antibody) at dilution of 1:50