

Anti-FUBP1/FBP Rabbit antibody



Catalog Number: 168920

Product name

Rabbit IgG

Anti-FUBP1/FBP Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

1:1000-1:5000

FUBP1/FBP Rabbit monoclonal antibody

1:50

Preparation

1:50


Antigen: A synthetic peptide of human FUBP1

1:20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Validations



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of FUBP1/FBP in Jurkat,C6 cell lysates using FUBP1/FBP Rabbit mAb(1:1000 diluted).Predicted band size:68kDa.Observed band size:74kDa.

Clonality

Monoclonal

Ig Type