

# NPR2 antibody



Catalog Number: 113220

## Product name

NPR2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

NPR2 Rabbit Polyclonal antibody. Positive WB detected in transfected cells. Observed molecular weight by Western-blot: 111 kDa, 117 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003995). Purification method: Antigen affinity purified.

## Formulation

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



WB result of anti-NPR2 (Catalog No:113220, 1:1000) with HEK293 overexpressed NPR2 by Jeremy Egbert and Laurinda Jaffe. (Another anti serum 6328 as control)