

# Anti-Liver FABP (8G6) Mouse antibody



Catalog Number: 168220

## Product name

Anti-Liver FABP (8G6) Mouse antibody

## Specificity

Human

## Antibody description

Liver FABP (8G6) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human L-FABP protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

## Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

ELISA

## Validations



Observed L-FABP levels in several Human urine samples and Human serum samples at different dilution.



Standard Curve for L-FABP: Capture Antibody Mouse mAb (8G6-D5) to L-FABP and Detector Antibody Mouse mAb(2D6-F6) to L-FABP.