

Human UBE2L6/UBCH8 (His tag) recombinant protein



Catalog Number: 519220

General Information

Gene Name Synonym

Ubiquitin/ISG15-Conjugating Enzyme E2 L6, Retinoic Acid-Induced Gene B Protein, RIG-B, UbcH8, Ubiquitin Carrier Protein L6, Ubiquitin-Protein Ligase L6, UBE2L6, UBCH8

Protein Construction

A DNA sequence encoding the Human UBE2L6/UBCH8 Met1-Ser153 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM HEPES, 100mM NaCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2L6/UBCH8 (His tag) recombinant protein