

SSX5 antibody



Catalog Number: 115620

Product name

SSX5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SSX5 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, A431 cells, mouse testis tissue. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of SSX5 recombinant protein (Accession Number: NM_021015). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115620(SSX5 antibody) at dilution of 1:300