

# P2RY11 Polyclonal Antibody



Catalog Number: 166620

---

**Product name**

Rabbit IgG

P2RY11 Polyclonal Antibody

**Applications****Specificity**

IHC

Human, Rat

**Dilutions****Antibody description**

IHC 1:50 - 1:100

Polyclonal antibody to P2RY11

**Validations****Preparation**

Antigen: Recombinant protein of human P2RY11

**Formulation**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

**Clonality**

Polyclonal

**Ig Type**

Immunohistochemistry - P2RY11 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human rectal cancer using P2RY11 antibody at dilution of 1:200 (40x lens).