

Anti-ARF5 Rabbit antibody



Catalog Number: 168620

Product name

Anti-ARF5 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ARF5 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human ARF5

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB: 1/1000

FC: 1/20

Validations



Western blot detection of ARF5 in Jurkat,C6 cell lysates using ARF5 Rabbit mAb(1:1000 diluted).Predicted band size:21 kDa.Observed band size:18kDa.