

ENOX1 antibody

Catalog Number: 110320

Product name

ENOX1 antibody

Specificity

Human; other species not tested.

Antibody description

ENOX1 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in Jurkat cells, MCF7 cells, Raji cells. Observed molecular weight by Western-blot: 73 kDa

Preparation

This antibody was obtained by immunization of ENOX1 recombinant protein (Accession Number: XM_024449374). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110320(ENOX1 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:110320(ENOX1 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:110320(ENOX1 Antibody) at dilution of 1:50