

GYPC antibody



Catalog Number: 111320

Product name

GYPC antibody

Specificity

Human; other species not tested.

Antibody description

GYPC Rabbit Polyclonal antibody. Positive WB detected in human blood tissue. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of GYPC recombinant protein (Accession Number: NM_002101). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:111320(GYPC Antibody) at dilution of 1:2000