

# Human ECH1 recombinant protein



Catalog Number: 518320

## General Information

### Gene Name Synonym

Delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase, mitochondrial,ECH1,, REF: C1008

### Protein Construction

A DNA sequence encoding the Human ECH1 Thr34-Leu328 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM Tris,100mM NaCl,10% Glycerol,pH 8.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human ECH1 recombinant protein