

Human IgA Heavy Chain antibody



Catalog Number: 111720

Product name

Human IgA Heavy Chain antibody

Specificity

Human; other species not tested.

Antibody description

Human IgA Heavy Chain Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue, human lung tissue. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of Human IgA Heavy Chain recombinant protein (Accession Number: BC016369). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

Validations



human plasma (0.3 ug) tissue were subjected to SDS PAGE followed by western blot with Catalog No:111720(Human IgA Antibody) at dilution of 1:1000. (1 Sec exposure)