

peroxiredoxin 2 antibody



Catalog Number: 113720

Product name

peroxiredoxin 2 antibody

Specificity

Human; other species not tested.

Antibody description

peroxiredoxin 2 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005809). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113720(PRD2 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using PRDX2 antibody Catalog No:113720 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).