

# RLBP1L2 antibody



Catalog Number: 114720

## Product name

RLBP1L2 antibody

## Immunogen

[Human RLBP1L2 Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

RLBP1L2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, Jurkat cells, L02 cells. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 38kd

## Preparation

This antibody was obtained by immunization of RLBP1L2 recombinant protein (Accession Number: NM\_001010852). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114720(RLBP1L2 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using RLBP1L2 antibody Catalog No:114720 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).