

STK33 antibody

Catalog Number: 115720

Product name

STK33 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

STK33 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, HeLa cells, mouse testis tissue. Positive IP detected in HEK-293 cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of STK33 recombinant protein (Accession Number: NM_001352388). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:115720(STK33 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:115720(STK33 antibody) at dilution of 1:100



+ - Normal Rabbit IgG
- + STK33 antibody

(under 10x lens)

IP Result of anti-STK33 (IP:Catalog No:115720,

STK33 antibody



Catalog Number: 115720

4ug; Detection:Catalog No:115720 1:300) with HEK-293 cells lysate 880ug.