

Anti-TNFRSF11B antibody



Catalog Number: 176720

Product name

Sheep IgG

Anti-TNFRSF11B antibody

Applications

Specificity

WB, IP, IF, IHC-P

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1:100-1:1,000

Sheep polyclonal antibody to TNFRSF11B

IP: 1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)

Preparation

IF: 1:50-1:500

This antigen of this antibody was peptide

IHC-P: 1:50-1:500

Formulation

Liquid, 1*TBS (pH7.4), 1%BSA, 40%Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C

Clonality

Polyclonal

Ig Type

Validations




Fig1: Western blot analysis of OPG monomer expression in rat bone marrow extract.