

# ACADM antibody



Catalog Number: 107720

## Product name

ACADM antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ACADM Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HEK-293 cells, HepG2 cells, K-562 cells, mouse heart tissue, rat heart tissue. Positive IHC detected in human liver tissue, human heart tissue. Observed molecular weight by Western-blot: 45kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000016). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107720(ACADM antibody) at dilution of 1:500