

Anti-Follistatin Rabbit antibody



Catalog Number: 168720

Product name

Anti-Follistatin Rabbit antibody

Specificity

Human, Rat

Antibody description

Follistatin Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Follistatin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

1:1000

Validations



Western blot detection of Follistatin in Rat Brain lysates using Follistatin Rabbit mAb(1:1000 diluted).Predicted band size:38kDa.Observed band size:38kDa.