

C3orf20 antibody

Catalog Number: 108720

Product name

C3orf20 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

C3orf20 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 101 kDa, 88 kDa

Preparation

This antibody was obtained by immunization of C3orf20 recombinant protein (Accession Number: NM_032137). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108720(C3orf20 antibody) at dilution of 1:500



IP Result of anti-C3orf20 (IP:Catalog No:108720, 3ug; Detection:Catalog No:108720 1:500) with mouse brain tissue lysate 10000ug.