

APOD antibody

Catalog Number: 108030

Product name

APOD antibody

Specificity

Human; other species not tested.

Antibody description

APOD Rabbit Polyclonal antibody. Positive WB detected in human serum tissue, HepG2 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 22-33 kDa

Preparation

This antibody was obtained by immunization of APOD recombinant protein (Accession Number: NM_001647). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human serum tissue were subjected to SDS PAGE followed by western blot with Catalog No:108030(APOD antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:108030(APOD antibody) at dilution of 1:100 (under 25x lens)