

# Anti-TGF beta 1 Rabbit antibody



Catalog Number: 169030

## Product name

Anti-TGF beta 1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

TGF beta 1 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of mouse TGF beta 1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

1:1000

## Validations



Western blot detection of TGF beta 1 in Jurkat,C6 cell lysates using TGF beta 1 Rabbit mAb(1:1000 diluted).Predicted band size:44kDa.Observed band size:44kDa.