

STRADA Polyclonal Antibody



Catalog Number: 160430

Product name

STRADA Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to STRADA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-210 of human STRADA (NP_699166.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:1000 - 1:2000

Validations



Western blot - STRADA Polyclonal Antibody

Western blot analysis of extracts of Mouse skeletal muscle, using STRADA antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.