

# MAPK10 antibody



Catalog Number: 112430

## Product name

MAPK10 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

MAPK10 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 45-48kd, 54-57kd

## Preparation

This antibody was obtained by immunization of MAPK10 recombinant protein (Accession Number: XM\_017008448). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112430(MAPK10 antibody) at dilution of



1:1000

Immunofluorescent analysis of HepG2 cells, using MAPK10 antibody Catalog No:112430 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).