

SMARCAL1 antibody



Catalog Number: 115430

Product name

SMARCAL1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SMARCAL1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Jurkat cells, mouse testis tissue. Observed molecular weight by Western-blot: 106kd

Preparation

This antibody was obtained by immunization of SMARCAL1 recombinant protein (Accession Number: XM_005246631). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115430(SMARCAL1 Antibody) at dilution of 1:600