

Anti-PTPRJ antibody

Catalog Number: 176430

Product name

Anti-PTPRJ antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to PTPRJ

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human ptpRJ. aa range:861-910

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1:: Western Blot analysis of HeLa, 293 cells using CD148 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).




Fig2:: Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:100