

# BCL7B antibody



Catalog Number: 108430

## Product name

BCL7B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BCL7B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue. Positive IHC detected in human colon cancer tissue, human lung cancer tissue, human pancreas tissue. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of BCL7B recombinant protein (Accession Number: NM\_001707). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117103(BCL7B antibody) at dilution of 1:400



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:117103(BCL7B Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:117103(BCL7B Antibody) at dilution of 1:50

# BCL7B antibody



Catalog Number: 108430

---

(under 40x lens)