

# CNNM4 antibody



Catalog Number: 109430

## Product name

CNNM4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CNNM4 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue, HEK-293 cells. Observed molecular weight by Western-blot: 90-100 kDa

## Preparation

This antibody was obtained by immunization of CNNM4 recombinant protein (Accession Number: NM\_020184). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109430(CNNM4 Antibody) at dilution of 1:600



IP Result of anti-CNNM4 (IP:Catalog No:109430, 4ug; Detection:Catalog No:109430 1:1000) with mouse testis tissue lysate 4400ug.