

# WIPI2 antibody



Catalog Number: 116830

## Product name

WIPI2 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

WIPI2 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, mouse testis tissue. Observed molecular weight by Western-blot: 49kd

## Preparation

This antibody was obtained by immunization of WIPI2 recombinant protein (Accession Number: NM\_015610). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:116830(WIPI2 antibody) at dilution of 1:200