

WIPI2 antibody



Catalog Number: 116830

Product name

WIPI2 antibody

Specificity

human, mouse; other species not tested.

Antibody description

WIPI2 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, mouse testis tissue. Observed molecular weight by Western-blot: 49kd

Preparation

This antibody was obtained by immunization of WIPI2 recombinant protein (Accession Number: NM_015610). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:116830(WIPI2 antibody) at dilution of 1:200