

# AGO2 antibody

Catalog Number: 107830

## Product name

AGO2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

AGO2 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in Jurkat cells, HEK-293 cells, K-562 cells. Observed molecular weight by Western-blot: 90-95kd

## Preparation

This antibody was obtained by immunization of AGO2 recombinant protein (Accession Number: BC007633). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107830(EIF2C2 antibody) at dilution of 1:300



IP Result of anti-AGO2 (IP:Catalog No:107830, 5ug; Detection:Catalog No:107830 1:500) with HEK-293 cells lysate 1800ug.