

# STK3 Polyclonal Antibody



Catalog Number: 165130

## Product name

STK3 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to STK3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 323-450 of human STK3 (NP\_006272.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

## Validations




Western blot - STK3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using STK3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Immunohistochemistry - STK3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using STK3 antibody at dilution of



1:100 (40x lens).


# STK3 Polyclonal Antibody

Catalog Number: 165130

---

Immunohistochemistry - STK3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using STK3 antibody at dilution of




1:100 (40x lens).

Immunohistochemistry - STK3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat

kidney using STK3 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - STK3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse stomach using STK3 antibody at dilution of 1:100 (40x lens).