

CD3 gamma antibody



Catalog Number: 107130

Product name

CD3 gamma antibody

Specificity

Human; other species not tested.

Antibody description

CD3 gamma Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 20-25 kDa

Preparation

This antibody was obtained by immunization of CD3 gamma recombinant protein (Accession Number: NM_000073). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:50-1:500


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107130(CD3G Antibody) at dilution of 1:4000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107130(CD3G Antibody) at dilution of 1:400 (under 10x lens). heat mediated antigen retrieved with Sodium Citrate buffer (pH6).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107130(CD3G Antibody) at dilution of 1:400 (under 40x lens). heat mediated antigen retrieved

CD3 gamma antibody



Catalog Number: 107130

with Sodium Citrate buffer (pH6).