

Mouse LY9/CD229/SLAMF3 (His Tag) recombinant protein



Catalog Number: 501530

General Information

Gene Name Synonym

Cell surface molecule Ly-9; Lymphocyte antigen 9; SLAM family member 3; Signaling lymphocytic activation molecule 3

Protein Construction

A DNA sequence encoding the extracellular domain of mouse Ly9 (NP_032560.2) (Met 1-Arg 453) was expressed, with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Lys 48

Molecular Mass

The recombinant mouse Ly9 consists of 417 amino acids and has a predicted molecular mass of 46.5 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rmLy9 is approximately 65 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse LY9 / CD229 / SLAMF3 Protein (His Tag)
SDS-PAGE