

H1N1 HA (His Tag) recombinant protein



Catalog Number: 502530

General Information

Protein Construction

A DNA sequence encoding the extracellular domain of influenza A virus hemagglutinin (A/Beijing/22808/2009 (H1N1)) (ADD64203.1) (Met 1-Gln 529) was fused with a polyhistidine tag at the C-terminus.

Organism

H1N1

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 18

Molecular Mass

The secreted recombinant hemagglutinin of H1N1(A/Beijing/22808/2009)HA (A/Beijing/22808/2009 (H1N1)) comprises 523

amino acids and has a predicted molecular mass of 59 kDa. As a result of glycosylation, it migrates as an approximately 72 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H1N1 (A/Beijing/22808/2009) Hemagglutinin / HA Protein (His Tag) SDS-PAGE