

PFKM Polyclonal Antibody



Catalog Number: 164530

Product name

PFKM Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to PFKM

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 501-780 of human PFKM (NP_001160160.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - PFKM Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PFKM antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.




Immunohistochemistry - PFKM Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using PFKM antibody at dilution of 1:100 (40x lens).

PFKM Polyclonal Antibody



Catalog Number: 164530



Immunohistochemistry - PFKM Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophageal cancer using PFKM antibody at dilution of 1:100 (40x lens).