

# Anti-Cytokeratin 9 Rabbit antibody



Catalog Number: 167530

## Product name

Anti-Cytokeratin 9 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Cytokeratin 9 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human Cytokeratin 9

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IP


## Dilutions

WB: 1/1000-1/2000

IHC: 1/50-1/200

IP: 1/50

## Validations



Western blot analysis of Cytokeratin 9 expression in HeLa cell lysate.