

MYH10-Specific antibody



Catalog Number: 112930

Product name

MYH10-Specific antibody

IP: 1:500-1:5000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MYH10-Specific Rabbit Polyclonal antibody. Positive FC detected in HepG2 cells. Positive IF detected in HepG2 cells. Positive WB detected in mouse brain tissue, human brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 229 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005964). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, FC, IP

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112930(MYH10-Specific antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells, using MYH10 antibody Catalog No:112930 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit




IgG (red).

MYH10-Specific antibody



Catalog Number: 112930

IP Result of anti-MYH10-Specific (IP:Catalog No:112930, 4ug; Detection:Catalog No:112930 1:1000) with mouse brain tissue lysate 4000ug.



1X10⁶ HepG2 cells were stained with 0.2ug MYH10-Specific antibody (Catalog No:112930, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.