

CD1B Polyclonal Antibody



Catalog Number: 164930

Product name

CD1B Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CD1B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 18-303 of human CD1B (NP_001755.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - CD1B Polyclonal Antibody

Western blot analysis of extracts of NCI-H125 cells, using CD1B Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.